



HAL
open science

Use of innovative technology for analysis of 54 VOCs in water: ITEX solution

Anne Togola, Benjamin Girardeau

► **To cite this version:**

Anne Togola, Benjamin Girardeau. Use of innovative technology for analysis of 54 VOCs in water: ITEX solution. SETAC Europe 21st Annual Meeting, May 2011, Milan, Italy. hal-00578154

HAL Id: hal-00578154

<https://brgm.hal.science/hal-00578154>

Submitted on 18 Mar 2011

HAL is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.

Use of innovative technology for analysis of 54 VOCs in water: ITEX solution

TOGOLA A., Girardeau B.

BRGM, 3, avenue Claude Guillemin, 45060 ORLEANS cedex 2, France,

E-mail contact: a.togola@brgm.fr

Abstract for a poster presentation

Topic : Environmental & analytical chemistry - B04 - Environmental fate and bioaccumulation of organic pollutants in aquatic systems

Key words : waste water, groundwater, VOCs analysis, degradation processes

Proposed chairs:

Contamination of aquifers by one or more volatile organic compounds (VOCs) is a national issue of potential concern because of the widespread and long-term use of many of these compounds. Many (VOCs) are toxic and some are considered to be carcinogenic, mutagenic, and teratogenic. Contamination also depends on the locations and types of VOC sources, the relative locations of wells, and the transport and fate of VOCs. The occurrence of low-level contamination of one or more VOCs in an aquifer also can provide managers with an early indication of the presence of VOCs that eventually might adversely affect the quality of water. Considering the high variability of VOCs concentration in groundwater especially in the case of industrial impacted groundwater, VOCs analysis is an important and sensitive point for understanding behaviour, attenuation or determining sources of the pollution.

Rapid sensitive and robust analysis of 54 VOCs including, BTEX and chlorinated solvents was developed using In Tube Extraction (ITEX) combined with GCMS for analysis in groundwater, effluent and surface waters. The study presents optimization of the method, validation on several matrices (occurrence of suspended, matter, organic matter content). Validations, considering potential trap needle saturation process and matrix effects were undertaken, highlighting robustness of the ITEX technology.

Detection limits ranged from 0.25 µg/L to 2 µg/L were validated with recoveries ranged from 85% to 103% depending on the studied compounds. Applications for the monitoring of VOCs in polluted areas will be presented showing the efficiency and accuracy of ITEX use and the obtaining of robust data concerning degradation of VOCs and natural attenuation for some of them.